

Product Datasheet

Nivolumab biosimilar antibody (orb1290228)

Description	Nivolumab biosimilar antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	PDCD1/PD1/CD279
Preservatives	PBS buffer pH 7.5
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG4-kappa
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	946414-94-4
Expiration Date	12 months from date of receipt.