

Product Datasheet

Rituximab biosimilar antibody (orb1290267)

Description	Rituximab biosimilar antibody
Species/Host	Mouse
Reactivity	Chimeric
Conjugation	Unconjugated
Tested Applications	FC
Target	MS4A1
Preservatives	PBS buffer pH 7.5
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1-kappa
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	174722-31-7
Expiration Date	12 months from date of receipt.