



Product Datasheet

4E-BP1 (phospho-Thr69) antibody (orb129028)

Description	Rabbit polyclonal antibody to 4E-BP1 (phospho-Thr69).
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13541, Q60876
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from J...



Immunohistochemistry analysis of formalin...