


Product Datasheet

Ras-GRF1 (phospho-Ser916) antibody (orb129030)

Description	Rabbit Polyclonal to Ras-GRF1 (phospho-Ser916).
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27671, Q13972
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...