



Product Datasheet

Caspase 9 (Cleaved-Asp353) antibody (orb129033)


Description	Rabbit polyclonal antibody to Caspase 9 (Cleaved-Asp353)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9R0S9, Q9JHK1
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of extracts from N...



Immunohistochemical analysis of formalin...