

Product Datasheet

Syndecan 2 peptide (orb12911)

Description	Syndecan 2 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Syndecan 2
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb13481 (Syndecan 2 antibody)
Expiration Date	6 months from date of receipt.