

Product Datasheet

FCRL5 Antibody (orb1291201)

Description	Anti-FCRL5 antibody(DMC220), IgG1 Chimeric mAb
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	FCRL5
Form/Appearance	Liquid
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	Rabbit/Human Fc chimeric IgG1
Clonality	Recombinant
Uniprot ID	Q96RD9
Dilution Range	ELISA 1/5000-10000 Flow Cyt 1/100
Expiration Date	12 months from date of receipt.