

Product Datasheet

M-CSF Antibody (orb1291216)

Description	Anti-M-CSF(Iacnotuzumab biosimilar) mAb
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	M-CSF
Form/Appearance	Powder
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	IgG1
Clonality	Recombinant
Uniprot ID	P09603
Dilution Range	ELISA 1:5000-10000 Flow Cyt 1:100
Expiration Date	12 months from date of receipt.