

Product Datasheet

Human UNC5B ELISA Kit (orb1291272)

Description	Human UNC5B (Netrin Receptor UNC5B) ELISA Kit
Reactivity	Human
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8IZJ1
Solubility (25°C)	10 ng/mL
Sensitivity	0.058 ng/mL
Expiration Date	Please enquire.

