

Product Datasheet

Rat FSHR ELISA Kit (orb1291276)

Description	Rat FSHR (Follicle Stimulating Hormone Receptor) ELISA Kit
Reactivity	Rat
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P20395
Solubility (25°C)	20 ng/mL
Sensitivity	0.124 ng/mL
Expiration Date	Please enquire.

