

Product Datasheet

Human CAF ELISA Kit (orb1291283)

Description	Human CAF (C-terminal fragement of agrin) ELISA Kit
Reactivity	Human
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	Q13111
Solubility (25°C)	2000 pg/mL
Sensitivity	13.1 pg/mL
Expiration Date	Please enquire.

