

Product Datasheet

Mouse dsDNA ELISA Kit (orb1291293)

Description	Mouse dsDNA (anti-double stranded DNA antibody IgM) ELISA KIT
Reactivity	Mouse
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, urine and other biological fluids
Solubility (25°C)	100 ng/mL
Sensitivity	0.58 ng/mL
Expiration Date	Please enquire.

