

Product Datasheet

Human TYRO3 ELISA Kit (orb1291297)

Description	Human TYRO3 (Tyrosine-protein kinase receptor TYRO3) ELISA Kit
Reactivity	Human
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids
Uniprot ID	Q06418
Solubility (25°C)	5000 pg/mL
Sensitivity	27 pg/mL
Expiration Date	Please enquire.

