

Product Datasheet

Rabbit ABCA1 ELISA Kit (orb1291300)

Description	Rabbit ABCA1 (ATP Binding Cassette Transporter A1) ELISA Kit
Reactivity	Rabbit
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, tissue homogenates and other biological fluids.
Solubility (25°C)	10 ng/mL
Sensitivity	0.06 ng/mL
Expiration Date	Please enquire.

