

Product Datasheet

Rabbit FASN ELISA Kit (orb1291303)

Description	Rabbit FASN (Fatty Acid Synthase) ELISA Kit
Reactivity	Rabbit
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P07855
Solubility (25°C)	10 ng/mL
Sensitivity	0.057 ng/mL
Expiration Date	Please enquire.

