

Product Datasheet

Rabbit CEBPa ELISA Kit (orb1291304)

Description	Rabbit CEBPa (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
Reactivity	Rabbit
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Tissue homogenates, cell lysates and other biological fluids
Solubility (25°C)	10 ng/mL
Sensitivity	0.065ng/mL
Expiration Date	Please enquire.

