

Product Datasheet

Rabbit Vaspin ELISA Kit (orb1291306)

Description	Rabbit Vaspin (Visceral Adipose Tissue Derived Serine Protease Inhibitor) ELISA Kit
Reactivity	Rabbit
Range	93.75-6000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma and other biological fluids.
Solubility (25°C)	6000 pg/mL
Sensitivity	35.19 pg/mL
Expiration Date	Please enquire.

