

Product Datasheet

Rabbit CPT1A ELISA Kit (orb1291309)

Description	Rabbit CPT1A (Carnitine Palmitoyltransferase 1A, Liver) ELISA Kit
Reactivity	Rabbit
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Tissue homogenates, cell lysates and other biological fluids
Solubility (25°C)	10 ng/mL
Sensitivity	0.062 ng/mL
Expiration Date	Please enquire.

