

Product Datasheet

Horse CD40L ELISA Kit (orb1291319)

Description	Horse CD40L (Cluster Of Differentiation 40 Ligand) ELISA Kit
Reactivity	Equine
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids
Solubility (25°C)	1000 pg/mL
Sensitivity	6.4 pg/mL
Expiration Date	Please enquire.

