

Product Datasheet

Horse CD8 ELISA Kit (orb1291320)

Description	Horse CD8 (Cluster of Differentiation 8) ELISA Kit
Reactivity	Equine
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Solubility (25°C)	100 ng/mL
Sensitivity	0.64 ng/mL
Expiration Date	Please enquire.

