

Product Datasheet

Sheep IL5 ELISA Kit (orb1291322)

Description	Sheep IL5 (Interleukin 5) ELISA Kit
Reactivity	Sheep
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q28586
Solubility (25°C)	1000 pg/mL
Sensitivity	5.7 pg/mL
Expiration Date	Please enquire.

