

Product Datasheet

Sheep IL13 ELISA Kit (orb1291323)

Description	Sheep IL13 (Interleukin 13) ELISA Kit
Reactivity	Sheep
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Solubility (25°C)	1000 pg/mL
Sensitivity	6.5 pg/mL
Expiration Date	Please enquire.

