

Product Datasheet

Human ACADL ELISA Kit (orb1291326)

Description	Human ACADL (Acyl Coenzyme A Dehydrogenase, Long Chain) ELISA Kit
Reactivity	Human
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P28330
Solubility (25°C)	10 ng/mL
Sensitivity	0.058 ng/mL
Expiration Date	Please enquire.

