

Product Datasheet

Rat FHL1 ELISA Kit (orb1291329)

Description	Rat FHL1 (Four And A Half LIM Domains Protein 1) ELISA Kit
Reactivity	Rat
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9WUH4
Solubility (25°C)	10 ng/mL
Sensitivity	0.058 ng/mL
Expiration Date	Please enquire.

