

Product Datasheet

Mouse ANA ELISA Kit (orb1291333)

Description	Mouse ANA (anti-nuclear Antibody) ELISA Kit
Reactivity	Human
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Solubility (25°C)	100 ng/mL
Sensitivity	0.58 ng/mL
Expiration Date	Please enquire.

