

Product Datasheet

Human KDM2B ELISA Kit (orb1291337)

Description	Human KDM2B (Lysine (K) specific demethylase 2B) ELISA kit
Reactivity	Human
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma and other biological fluids
Uniprot ID	Q8NHM5
Solubility (25°C)	20 ng/mL
Sensitivity	0.115 ng/mL
Expiration Date	Please enquire.

