

Product Datasheet

Human P4HTM ELISA Kit (orb1291359)

Description	Human P4HTM (Transmembrane prolyl 4-hydroxylase) ELISA Kit
Reactivity	Human
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9NXG6
Solubility (25°C)	1000 pg/mL
Sensitivity	5.7 pg/mL
Expiration Date	Please enquire.

