

## Product Datasheet

### Human TNFRSF8 ELISA Kit (orb1291360)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Human TNFRSF8 (Tumor Necrosis Factor Receptor Superfamily, Member 8) ELISA Kit |
| <b>Reactivity</b>          | Human                                                                          |
| <b>Range</b>               | 0.32-20 ng/mL                                                                  |
| <b>Tested Applications</b> | ELISA                                                                          |
| <b>Note</b>                | For research use only                                                          |
| <b>Application notes</b>   | Sandwich                                                                       |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates and other biological fluids                  |
| <b>Uniprot ID</b>          | <b>P28908</b>                                                                  |
| <b>Solubility (25°C)</b>   | 20 ng/mL                                                                       |
| <b>Sensitivity</b>         | 0.135 ng/mL                                                                    |
| <b>Expiration Date</b>     | Please enquire.                                                                |

