

Product Datasheet

Rat nAChR ELISA Kit (orb1291366)

Description	Rat nAChR (Nicotinic Acetylcholine Receptor) ELISA Kit
Reactivity	Rat
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, tissue homogenates and other biological fluids.
Solubility (25°C)	10 ng/mL
Sensitivity	0.059 ng/mL
Expiration Date	Please enquire.

