

Product Datasheet

Rat VEGFD ELISA Kit (orb1291367)

Description Rat VEGFD (Vascular Endothelial Growth Factor D) ELISA Kit

Reactivity Rat

Range 15.63-1000 pg/mL

Tested Applications ELISA

Note For research use only

Application notes Sandwich

Sample Types serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Uniprot ID [O35251](#)

Solubility (25°C) 1000 pg/mL

Sensitivity 5.7 pg/mL

Expiration Date Please enquire.

