

Product Datasheet

Rat AANAT ELISA Kit (orb1291368)

Description Rat AANAT (Arylalkylamine-N-Acetyltransferase) ELISA Kit

Reactivity Rat

Range 0.16-10 ng/mL

Tested Applications ELISA

Note For research use only

Application notes Sandwich

Sample Types tissue homogenates, cell lysates and other biological fluids

Uniprot ID [Q64666](#)

Solubility (25°C) 10 ng/mL

Sensitivity 0.059 ng/mL

Expiration Date Please enquire.

