

Product Datasheet

Rabbit BALP ELISA Kit (orb1291369)

Description	Rabbit BALP (BoneAlkaline Phosphatase) ELISA Kit
Reactivity	Rabbit
Range	0.79-50 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates and other biological fluids
Solubility (25°C)	50 ng/mL
Sensitivity	0.3 ng/mL
Expiration Date	Please enquire.

