

Product Datasheet

Mouse TNFSF12 ELISA Kit (orb1291372)

Description	Mouse TNFSF12 (Tumor Necrosis Factor Ligand Superfamily, Member 12) ELISA Kit
Reactivity	Mouse
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	O54907
Solubility (25°C)	1000 pg/mL
Sensitivity	6.9 pg/mL
Expiration Date	Please enquire.

