

Product Datasheet

Human YWHAG ELISA Kit (orb1291376)

Description Human YWHAG (14-3-3 protein gamma) ELISA Kit

Reactivity Human

Range 0.32-20 ng/mL

Tested Applications ELISA

Note For research use only

Application notes Sandwich

Sample Types Serum, plasma, tissue homogenates and other biological fluids.

Uniprot ID **P61981**

Solubility (25°C) 20 ng/mL

Sensitivity 0.131 ng/mL

Expiration Date Please enquire.

