

Product Datasheet

Human F-PSA ELISA Kit (orb1291378)

Description	Human F-PSA (Free ProstateSpecific Antigen) ELISA Kit
Reactivity	Human
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Solubility (25°C)	2000 pg/mL
Sensitivity	17.75 pg/mL
Expiration Date	Please enquire.

