

Product Datasheet

Human TRIM72 ELISA Kit (orb1291379)

Description	Human TRIM72 (Tripartite motif-containing protein 72) ELISA kit
Reactivity	Human
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q6ZMU5
Solubility (25°C)	10 ng/mL
Sensitivity	0.097 ng/mL
Expiration Date	Please enquire.

