

Product Datasheet

NF-L antibody (orb1291513)

Description	NF-L antibody
Species/Host	Mouse
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human NF-L
Target	NF-L
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	68 kDa
Uniprot ID	P07196
Dilution Range	IHC-P=1:200-800
Expiration Date	12 months from date of receipt.