



Product Datasheet

Tau-4 antibody (orb1291514)


Description	Tau-4 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Human Tau-4 Protein (1-441/441 aa)
Target	Tau-4
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	46/79 kDa
Uniprot ID	P10636
Dilution Range	WB=1:500-2000, IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P analysis of rat cerebellum using T...



IHC-P analysis of mouse cerebellum using...



Western blot: Lane 1: Mouse Cerebellum; Lan...