



Product Datasheet

MRE11 antibody (orb1291515)


Description	MRE11 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human MRE11
Target	MRE11
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	80 kDa
Uniprot ID	P49959
Dilution Range	WB=1:500-1000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse testis
using Mre...



IHC-P
analysis of
rat spleen
using
Mre11...



IHC-P
analysis of
human tonsil
using Mre...