

Product Datasheet

Tau (phospho-Thr231) antibody [Out of stock] (orb1291593)

Description	Tau (phospho-Thr231) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr231 (VR(p-T)PP)
Target	phospho-Tau (Thr231)
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46/79 kDa
Uniprot ID	P10636
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.