




Product Datasheet

TRMT61B Antibody (orb1291636)

Description	TRMT61B Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, FC, IF, IHC, WB	
Immunogen	E.coli-derived human TRMT61B recombinant protein (Position: R7-N477).	IHC analysis of TRMT61B using anti-TRMT61B...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	IHC analysis of TRMT61B using anti-TRMT61B...
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	57 kDa	
Uniprot ID	Q9BVS5	
Dilution Range	Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10 ⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human	IHC analysis of TRMT61B using anti-TRMT61B...
Expiration Date	12 months from date of receipt.	