





Product Datasheet

TRMT6 Antibody (orb1291639)

Description	TRMT6 Antibody	
Species/Host	Rabbit	  IHC analysis of TRMT6 using anti-TRMT6 a...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, FC, ICC, IF, IHC, WB	
Immunogen	E.coli-derived human TRMT6 recombinant protein (Position: E6-K489).	 IHC analysis of TRMT6 using anti-TRMT6 a...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.	
Note	For research use only	 IHC analysis of TRMT6 using anti-TRMT6 a...
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	60 kDa	 IHC analysis of TRMT6 using anti-TRMT6 a...
Uniprot ID	Q9UJA5	
Dilution Range	Western blot, 0.1-0.25 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 1-2 µg/ml, Human Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10 ⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human	
Expiration Date	12 months from date of receipt.	