

Product Datasheet

Human CLDN18.2 full length protein-VLP (orb1291767)

Description	Human Full length Claudin 18.2 Protein-VLP
Reactivity	Human
Target	CLDN18.2
Form/Appearance	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at - 80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Source	Mammalian
MW	The Human Claudin 18.2 Protein has a predicted MW of 29 kDa. Due to PTM, the actual MW on SDS-PAGE gel is around 50kDa.
Uniprot ID	P56856
Expiration Date	6 months from date of receipt.