

Product Datasheet

Human CCR8 full length protein-MNP (orb1291769)

Description	Human CCR8 full length protein membrane nanoparticles (MNPs)
Reactivity	Human
Target	CCR8
Form/Appearance	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Source	Mammalian
MW	The human full length CCR8 Protein has a MW of 40.7 kDa
Uniprot ID	P51685
Expiration Date	6 months from date of receipt.