

Product Datasheet

Biotinylated Human IL-17 RC Protein (orb1291791)

Description Biotinylated Human IL-17 RC, Fc, is expressed from human 293 cells (HEK293). It contains AA Leu 21 -...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Biotin

Target IL-17 RC

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag C-Fc & C-Avi

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 77.6 kDa. The protein migrates as 80-110 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity 95%


Source Biotinylated Human IL-17 RC, Fc, (orb1291791) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - His 465 (Accession # Q8NAC3-2).

MW 77.6 kDa

Uniprot ID [Q8NAC3](#)

NCBI [NP_703190.2](#)

Expiration Date 6 months from date of receipt.



Immobilized
Human
IL17A, Tag
Free at 5
◆...