

Product Datasheet

CYP11B1 peptide (orb12918)

Description	CYP11B1 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	CYP11B1
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb5937 (CYP11B1 antibody)
Expiration Date	6 months from date of receipt.