

Product Datasheet

SARS-CoV-2 Nucleoprotein [Omicron B.1.1.529] protein [Out of stock] (orb1292412)

Description	Recombinant SARS-CoV-2 Nucleoprotein [Omicron B.1.1.529] protein
Reactivity	Virus
Tested Applications	ELISA, WB
Preservatives	Lyophilized from a 0.2 µm filtered solution in PBS with 5% trehalose, pH7.4
Storage	-20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution
Tag	N-terminal His Tag
Note	For research use only
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Formula	Lyophilized from a 0.2 µm filtered solution in PBS with 5% trehalose, pH _{7.4}
Purity	Greater than 95% by SDS-PAGE gel analyses
Source	E.Coli
MW	45.6 kDa
Uniprot ID	P0DTC9
Expiration Date	6 months from date of receipt.