

Product Datasheet

Substance P antibody antibody (orb1292438)

Description	Rabbit Polyclonal antibody to Substance P
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	KLH conjugated synthetic peptide derived from human Substance P
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	P20366
Dilution Range	IHC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.