

Product Datasheet

LAMP1 antibody (orb1292487)

Description	Mouse Monoclonal antibody to LAMP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	lysosomal-associated membrane protein 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	120 kDa
Uniprot ID	P11279
Dilution Range	WB: 1:1000 - 1:10000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.