


Product Datasheet

Rabbit CPT1B ELISA Kit (orb1292599)

Description	Rabbit CPT1B (Carnitine palmitoyltransferase 1B) ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/ml
Target	CPT1B
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	G1TW00
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rabbit
CPT1B(Carnitine
palmitoyltransfer...